

MICROTITER TEMPERATURE CONTROL PLATE

Uniform Surface Temperature at an Unbeatable Price



SSCS Part Number 14-13129-2

- Accepts all ANSI SBS 2004-1 microwell plates
- Designed for high heat transfer and unparalleled temperature uniformity of $<0.1^{\circ}$ C across plate surface
- Ready for use between -5° and 50° C
- Precision construction means no leveling required for fast and accurate installation
- Robotic alignment marks standard
- Insulated design minimizes unwanted condensation
- Comes with no drip twist lock fittings
- Insulated hose kits also available
- RoHS Compliant
- Value priced

Solid State
COOLING SYSTEMS
Temperature Control...Precisely.