

Recirculating Bath

Low cost bath accessory which works with any of our thermoelectric recirculating chillers



Recirculating Bath Accessory
Part Number: 10-14920-1
Serial Number: 9800168R
This bath accessory should only be used with
Solid State Cooling Systems chillers
Solid State Cooling Systems
Wappingers Falls, NY 12590
(845) 296-1300 www.sscooling.com

Turn any recirculating chiller into a temperature controlled bath

- Holds precise temperature to $\pm 0.1^{\circ}\text{C}$
- Tank volume is 10.6 liters (2.8 gallons)
- Easy to set up, easy to use

The Right Temperature Without Fail



Open bath with removable cover

Low cost recirculating bath

Instead of purchasing a standalone recirculating bath chiller, this simple bath accessory allows any existing recirculating chiller to be used for open bath applications. Compatible with any of our thermoelectric recirculating chillers including ThermoCube, ThermoCube II, TCube edge, and ThermoRack. Available with fittings on long side or short side. Three-layer insulated hoses sold separately.

SPECIFICATIONS

Operating Range	-10°C to 65°C
Ambient Temperature	10°C to 40°C non-condensing
Stability / Repeatability	±0.1°C at constant load
Bath Capacity	10.6 liters (2.8 US gallons)
Fluid Fittings	1/4" valved female CPC
Wetted Materials	Stainless steel and polymers
Size (L x W x H)	12.25" x 10" x 10" (31.1 x 25.4 x 25.4 cm)
Weight (empty)	5 lbs (2.4 kg)
Input Power	No power required
Warranty	2 years

Solid State Cooling Systems